



PLASTIQLINE
an All Care brand

Waste bin 20 litre with knee opening plastic, PQQBL

In an environment where hygiene is crucial, the PlastiQline plastic waste bin earns pride of place. The ingenious knee operation make the transfer of bacteria via the hands a thing of the past. In hospitals and care institutions hygiene is of utmost importance. Any transfer of bacteria must be prevented at all cost. The knee-operated PlastiQline plastic waste bin prevents any undesired unhygienic contamination. How does it work? The minute you press your knee against the lower container, the lid will pop open. Once the waste is deposited and the pressure on the lower container removed, the lid will pop back in its original position automatically. Do you attach great value to optimal hygiene? Then request a quote for the PlastiQline Knee-operated 20 Litre Plastic Waste Bin.

Technical description

Waste bin for wall mounting, knee opening system. The lid opens, when the knee touched the bin. With internal removable bin, with flats for fixing the waste bag.

Specificaties

Article number	5666
Material	Plastic (White)
H x W x D	520 mm x 350 mm x 250 mm
Weight	3 kg
Guarantee	1 year from the date of purchase
Application	Knee opening
Content	20 litre

All Care, manufacturer of dispenserlines