



Waste bin 8 litre with knee opening plastic, PQKBS

Hygiene plays an essential role in hospitals and care institutions where the transfer of bacteria must be prevented at all costs. When using the PlastiQline Plastic Waste Bin - fully knee operated - you will prevent any unwanted unhygienic contamination. It's very simple; the lid opens as soon as you press your knee against the bin. Once you have put the waste in the bin and stop pressing against it, the lid will automatically close. How valuable is optimal hygiene to you? Request a quote for the PlastiQline Plastic Waste Bin!

Technical description

Waste bin for wall mounting, knee opening system. The lid opens, when the knee touched the bin.With internal removable bin, with flats for fixing the waste bag.

Specificaties

Article number	5665
Material	Plastic (White)
H x W x D	360 mm x 220 mm x 220 mm
Weight	1,5 kg
Guarantee	1 year from the date of purchase
Application	Knee opening
Content	8 litre

All Care, manufacturer of dispenserlines

All Care B.V. Storkstraat 41 3905 KX Veenendaal P.O. Box 374 3900 AJ Veenendaal The Netherlands T : +31 (0) 318 517 778 F : +31 (0) 318 514 498 E : info@all-care.eu

